

# Multi-Drug Urine Test Cup

## INSTRUCTIONS FOR USE

PLEASE READ ALL INFORMATION IN THE INSTRUCTIONS FOR USE BEFORE USING THE TEST!

REF See Box Label



**Ketamine (KET)**  
Ketamine was developed in the 1960s to replace phenylcycidine (PCP) as an anesthetic agent and is most commonly used in veterinary medicine today. In addition to rohypnol (add hyperlink to page) and GHB, it is also considered a club drug and may be used in drug-facilitated sexual assault situations. It is odorless, tasteless and usually swallowed in powder form or injected. Once taken, it is very short-acting and shows effects within minutes. Under federal law, ketamine is classified as a Schedule III drug, meaning it has approved medical use, but still possesses a high potential for abuse.

**Methamphetamine (MET)**  
Methamphetamine is a potent sympathomimetic agent with therapeutic applications. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and a sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behavior, and cardiac dysrhythmias. The pattern of psychosis which may appear at half-life of about 15 hours is excreted in urine as amphetamine and oxidized as deaminated and hydroxylated derivatives. However, 40% of methamphetamine is excreted unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

**Methadone (MTD)**  
Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiate dependence (Heroin, Vicodin, Percocet, Morphine). It is administered either orally, or by intravenous or intra-muscular injection. The duration of effect of methadone is 12–24 hours. Its major urinary excretion products are methadone, EDDP (2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine), and EMDP (2-ethyl-5-methyl-3, 3-diphenylpyrrolidine).

**Methadone Metabolite (EDDP)**  
EDDP(2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine) is the primary metabolite of methadone. Methadone is a synthetic analgesic drug that is originally used in the treatment of narcotic addicts. The detection of EDDP is more beneficial than traditional methadone screening since EDDP exists only in urine from individuals that ingested methadone. The tampering of specimens by spiking the urine with methadone can be prevented. Secondly, renal clearance of EDDP is not affected by urinary pH, therefore the EDDP test provides a more accurate result of methadone ingestion than the methadone parent screen.

**Methylenedioxyamphetamine - ecstasy (MDMA)**  
MDMA belongs to a family of man-made drugs. Its relatives include MDA (methylenedioxyamphetamine), and MDEA (methylenedioxyethylamphetamine). They all share the amphetamine-like effects. MDMA is a stimulant with hallucinogenic tendencies described as an empathogen as it releases mood-altering chemicals, such as cartoning and L-Dopa, and may generate feelings of love and friendliness. The adverse effects of MDMA use include elevated blood pressure, hyperthermia, anxiety, paranoia and insomnia. MDMA is administered either by oral ingestion or intravenous injection. The effects of MDMA begin 30 minutes after intake, peak in an hour and last for 2–3 hours.

**Morphine (MOP)**  
Opiate refers to any drug that is derived from the opium poppy, including the natural products, morphine and codeine, and the semi-synthetic drugs such as heroin. Opioid is more general, referring to any drug that acts on the opioid receptor. Opioid analgesics comprise a large group of substances which control pain by depressing the central

This package insert applies to any combination of multi-drug tests and adulteration control tests. Therefore, some information on the performance characteristics of the product may not be relevant to your test. Please refer to the labels on the packaging and the prints on the test device to identify which drugs and adulteration controls are included in your test.

### INTENDED USE

Multi-Drug Urine Test Cup is a rapid urine screening test. It's a lateral flow, one-step immunoassay for the qualitative detection of specific drugs and their principal metabolites in human urine at specified cut-off concentrations, with additional semi quantitative adulteration controls. The multi-drug test device can be combined with the adulteration controls such as Creatinine (CRE), Glutaraldehyde (GLU), Nitrite (NIT), pH, Specific Gravity (S.G.), and/or Oxidants/Pyridinium Chlorochromate (OXI/PCC), which is used for the determination of diluted or adulterated urine specimens. The adulteration control is an important pre-screening test for drug-testing.

Drug (Identifier)	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP 1000)	d-Amphetamine	1000
Amphetamine (AMP 500)	d-Amphetamine	500
Amphetamine (AMP 300)	d-Amphetamine	300
Barbiturates (BAR 300)	Secobarbital	300
Benzodiazepines (BZO 300)	Oxazepam	300
Benzodiazepines (BZO 200)	Oxazepam	200
Benzodiazepines (BZO 100)	Oxazepam	100
Buprenorphine (BUP 10)	Buprenorphine	10
Cocaine (COC 300)	Benzoylcocaine	300
Cocaine (COC 150)	Benzoylcocaine	150
Cocaine (COC 100)	Benzoylcocaine	100
Cannabinoids (THC 50)	11-nor-Δ <sup>9</sup> -THC-9-COOH	50
Cannabinoids (THC 40)	11-nor-Δ <sup>9</sup> -THC-9-COOH	40
Cannabinoids (THC 25)	11-nor-Δ <sup>9</sup> -THC-9-COOH	25
Cotinine (COT)	Cotinine	200
Ethyl Glucuronide (ETG)	Ethyl Glucuronide	500
Fentanyl (FTY)	Fentanyl	20
Ketamine (KET 1000)	Ketamine	1000
Ketamine (KET 500)	Ketamine	500
Methamphetamine (MET 1000)	d-Methamphetamine	1000
Methamphetamine (MET 500)	d-Methamphetamine	500

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

Methamphetamine (MET 300)	d-Methamphetamine	300
Methadone (MTD 300)	Methadone	300
Methadone (MTD 200)	Methadone	200
Methadone Metabolite (EDDP 300)	2-ethylidene-1,5-dimethyl-3, 3-diphenylpyrrolidine (EDDP)	300
Methylenedioxyamphetamine -ecstasy (MDMA 500)	3,4-Methylenedioxyamphet-amine (MDMA)	500
Methylenedioxyamphetamine -ecstasy (MDMA 300)	3,4-Methylenedioxyamphet-amine HCl (MDMA)	300
Morphine (MOP 300)	Morphine	300
Morphine (MOP 100)	Morphine	100
Opiate (OPI)	Morphine	2000
Oxycodone (OXY)	Oxycodone	100
Phencyclidine (PCP)	Phencyclidine	25
Propoxyphene (PPX)	d-Propoxyphene	300
Synthetic Cannabis (K2 50)	JWH-018 / JWH-073	50
Synthetic Cannabis (K2 25)	JWH-018 / JWH-073	25
Tricyclic Antidepressants (TCA)	Nortriptyline	1000
Tramadol (TRA 1000)	Tramadol	1000
Tramadol (TRA 200)	Tramadol	200
6-Monoacetylmorphine (6-MAM)	6-Monoacetylmorphine	10

Conjugations of the Multi-Drug Urine Test Cup can consist of any combination of the above listed drug analytes. It is intended for forensic use only. This assay provides a qualitative, preliminary test result. A more specific analytical method must be used in order to obtain a confirmed result. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS) are the preferred confirmatory methods. Professional judgment should be applied to any drug test result, particularly when preliminary positive results are indicated.

### SUMMARY

**Amphetamine (AMP)**  
Amphetamines are a structurally related "designer" drugs are sympathomimetic amines whose biological effects include potent central nervous system (CNS) stimulation, anorectic, hyperthermic, and cardiovascular properties. They are usually taken orally, intravenously, or by smoking. Amphetamines are readily absorbed from the gastrointestinal tract and are then either deactivated by the liver or excreted unchanged in the urine. Methamphetamine is partially metabolized to amphetamine and its major active metabolite. Amphetamines increase the heart rate and blood pressure, and suppress the appetite. Some studies indicate that heavy abuse may result in permanent

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

damage to certain essential nerve structural in the brain. The effects of Amphetamines generally last 2–4 hours following use and the drug has a half-life of 4–24 hours in the body. About 30% of amphetamines are excreted in the urine in unchanged form, with the remainder as hydroxylated and deaminated derivatives. It can be detected in the urine for 1 to 2 days after use.

**Barbiturates (BAR)**  
Barbiturates are central nervous system depressants. They are usually administered orally but are sometimes injected intramuscularly and intravenously. Barbiturates range from short-acting (approximately 15 minutes, such as secobarbital) to long-acting (24 hours or longer, such as Phenobarbital). Short-acting barbiturates are extensively metabolized in the body, while the long-acting ones are secreted primarily unchanged. Barbiturates produce alertness, wakefulness, increased energy, reduced hunger, and an overall feeling of well being. Large doses of Barbiturate could develop tolerance and physiological dependency and lead to its abuse.

**Benzodiazepines (BZO)**  
Benzodiazepines are a class of drugs that are often therapeutically used as anxiolytics, anti-convulsants and sedative hypnotics. Benzodiazepines manifest their presence by analgesia, drowsiness, confusion, diminished reflexes, lowering of body temperature, respiratory depression, blockade of adrenergic response, and a decrease in peripheral resistance without an impact on the cardiac index. The major pathways of elimination are the kidneys (urine) and the liver where it is conjugated to glucuronic acid. Large doses of Benzodiazepines could develop tolerances and physiological dependency and lead to its abuse. Only trace amounts (less than 1%) of Benzodiazepines are excreted unaltered in the urine, most of Benzodiazepines in urine is conjugated drug. Oxazepam, a common metabolite of many benzodiazepines, remains detectable in urine for up to one week, which makes Oxazepam a useful marker of Benzodiazepines abuse.

**Buprenorphine (BUP)**  
Buprenorphine is a potent analgesic often used in the treatment of opioid addiction. The drug is sold under the trade names Subutex®, Buprenex®, Temgesic®, and Suboxone®, which contain Buprenorphine HCl alone or in combination with Naloxone HCl. Therapeutically, Buprenorphine is used as a substitution treatment for opioid addicts. Substitution treatment is a form of medical care offered to opiate addicts (primarily heroin addicts) based on a similar or identical substance to the drug normally used. In substitution therapy, Buprenorphine is as effective as Methadone but demonstrates a lower level of physical dependence. Concentrations of free Buprenorphine and Nortbuprenorphine in urine may be less than 1 ng/ml after therapeutic administration, but can range up to 20 ng/ml in abuse situations. The plasma half-life of Buprenorphine is 2–4 hours. While complete elimination of a single dose of the drug can take as long as 6 days, the window of detection for the parent drug in urine is thought to be approximately 3 days. Substantial abuse of Buprenorphine has also been reported in many countries where various forms of the drug are available. The drug has been diverted from legitimate channels through theft, doctor shopping, and fraudulent prescriptions, and been abused via intravenous, sublingual, intranasal and inhalation routes.

**Cocaine (COC)**  
Cocaine derived from leaves of coca plant, is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using cocaine are euphoria, confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Cocaine is excreted in urine primarily as benzoylecgonine in a short period of time.

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

**Cannabinoids (THC)**  
Cannabinoids is a hallucinogenic agent derived from the flowering portion of the hemp plant. The active ingredients in Cannabinoids, THC & Cannabinol can be metabolized and excreted as 11-nor-Δ<sup>9</sup>-tetrahydrocannabinol-9-carboxylic acid with a half-life of 24 hours. It can be detected for 1 to 5 days after use. Smoking is the primary method of use of Cannabinoids/cannabis. Higher doses used by abusers produce central nervous system effects, altered mood and sensory perceptions, loss of coordination, impaired short-term memory, anxiety, paranoia, depression, confusion, hallucinations and increased heart rate. A tolerance to the cardiac and psychotropic effects can occur, and withdrawal syndrome produces restlessness, insomnia, anxiety and nausea.

**Cotinine (COT)**  
Cotinine is the first-stage metabolite of nicotine, a toxic alkaloid that produces stimulation of the autonomic ganglia and central nervous system in humans. Nicotine is a drug to which virtually every member of a tobacco-smoking society is exposed whether through direct contact or second-hand inhalation. In addition to tobacco, nicotine is also commercially available as the active ingredient in smoking replacement therapies such as nicotine gum, transdermal patches and nasal sprays. In a 24-hour urine, approximately 5% of a nicotine dose is excreted as unchanged drug with 10% as cotinine and 85% as hydroxycotinine. The concentrations of other metabolites are believed to account for less than 5%. While cotinine is thought to be an inactive metabolite, its elimination profile is more stable than that of nicotine which is largely urine pH dependent. As a result, cotinine is considered a good biological marker for determining nicotine use. The plasma half-life of nicotine is approximately 60 minutes following inhalation or parenteral administration. Nicotine and cotinine are rapidly eliminated by the kidney; the window of detection for cotinine in urine at a cutoff level of 200 ng/mL is expected to be up to 2–3 days after nicotine use.

**Ethyl Glucuronide (ETG)**  
Ethyl Glucuronide (ETG) is a direct metabolite of alcohol. Presence in urine may be used to detect recent alcohol intake, even after alcohol is no longer measurable. Traditional laboratory methods detect the actual alcohol in the body, which reflects current intake within the past few hours (depending on how much was consumed). The presence of ETG in urine is a definitive indicator that it can be detected in the urine for 3 to 4 days after drinking alcohol, even alcohol is eliminated from the body. Therefore, ETG is a more accurate indicator of the recent intake of alcohol than measuring for the presence of alcohol itself. The ETG test can aid in the diagnosis of drunk driving and alcoholism, which has implications for the forensic identification and medical examination.

**Fentanyl (FTY)**  
Fentanyl belongs to powerful narcotics analgesics, and is a μ opioid receptors receptor stimulant. Fentanyl is one of the varieties that been listed in management of United Nations "Single Convention of narcotic drug in 1961". Among the opiates agents that under international control, fentanyl is one of the most commonly used to cure moderate to severe pain1. After continuous injection of fentanyl, the sufferer will have the performance of protracted opioid abstinence syndrome, such as ataxia and irritability etc2,3, which presents the addiction after taking fentanyl in a long time. Compared with drug addicts of amphetamine, drug addicts who take fentanyl mainly have got the possibility of high infective rate of HIV, more dangerous injection behavior and more lifelong medication overdose.

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19